

## DISPOSABLE STOPCOCK

R<sub>X</sub> Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

### Intended Use/Indications:

The Stopcock may be used for infusion of fluids and/or obtaining blood samples or other fluids.

### Contraindications:

The Stopcock is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for the cardiopulmonary bypass in general.



### Precautions:

**Caution:** Inspect carefully for leaks before and during use.

**Caution:** A leak at the connection of the Stopcock may cause fluid or blood loss, air emboli, or a loss of sterility.

**Caution:** The Stopcock is disposable. Do not reuse, reprocess or re-sterilize.

### Direction for Use:

1. Aseptically remove the Stopcock from the pouch and transfer it to the sterile field.
2. Connect the Stopcock to the desired component using a clockwise push-twist motion to secure into the position.
3. The Stopcock handle should be in position over the port to be closed. The remaining ports will be open.
4. Make certain the connection is secure and free of air.



**NovoSci™**

2021 Airport Rd.  
Conroe, TX 77301  
Toll Free 1-800-854-0567  
Direct 1-281-363-4949

424230-00A Rev. C

## DISPOSABLE STOPCOCK

R<sub>X</sub> Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

### Intended Use/Indications:

The Stopcock may be used for infusion of fluids and/or obtaining blood samples or other fluids.

### Contraindications:

The Stopcock is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for the cardiopulmonary bypass in general.



### Precautions:

**Caution:** Inspect carefully for leaks before and during use.

**Caution:** A leak at the connection of the Stopcock may cause fluid or blood loss, air emboli, or a loss of sterility.

**Caution:** The Stopcock is disposable. Do not reuse, reprocess or re-sterilize.

### Direction for Use:

1. Aseptically remove the Stopcock from the pouch and transfer it to the sterile field.
2. Connect the Stopcock to the desired component using a clockwise push-twist motion to secure into the position.
3. The Stopcock handle should be in position over the port to be closed. The remaining ports will be open.
4. Make certain the connection is secure and free of air.



**NovoSci™**

2021 Airport Rd.  
Conroe, TX 77301  
Toll Free 1-800-854-0567  
Direct 1-281-363-4949

424230-00A Rev. C