

DISPOSABLE PERFUSION ADAPTER CONNECTOR

R_x Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Intended Use/Indications:

The Perfusion Adapter Connectors are indicated whenever it is advisable to terminate perfusion tubing (1/4" I.D. or 3/8" I.D.) with a male luer taper.

Contraindications:

The Perfusion Adapter Connectors is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for perfusion in general.



Cautions and Warnings:

Caution: Insert the connector into the tubing over the third barn for security.

Caution: Check luer connection periodically for security.

Caution: The connector is disposable. Do not reuse or re-sterilize.

Caution: Contents of unopened, undamaged package are sterile.

Caution: Inspect packaging for damage prior to use. Sterility and/or product performance may be compromised if damage is found.

Direction for Use:

1. Aseptically remove the Adapter Connector from the pouch.
2. Insert the barbed end of the connector securely into the perfusion tubing.
3. Insert the male luer end of the connector into the desired female luer with a slight twisting motion to insure a secure tight seal.
4. De-bubble the line if it is intended for perfusion.



NovoSci[™]

2021 Airport Rd.
Conroe, TX 77301
Toll Free 1-800-854-0567
Direct 1-281-363-4949

424006-000 Rev. D

DISPOSABLE PERFUSION ADAPTER CONNECTOR

R_x Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Intended Use/Indications:

The Perfusion Adapter Connectors are indicated whenever it is advisable to terminate perfusion tubing (1/4" I.D. or 3/8" I.D.) with a male luer taper.

Contraindications:

The Perfusion Adapter Connectors is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for perfusion in general.



Cautions and Warnings:

Caution: Insert the connector into the tubing over the third barn for security.

Caution: Check luer connection periodically for security.

Caution: The connector is disposable. Do not reuse or re-sterilize.

Caution: Contents of unopened, undamaged package are sterile.

Caution: Inspect packaging for damage prior to use. Sterility and/or product performance may be compromised if damage is found.

Direction for Use:

1. Aseptically remove the Adapter Connector from the pouch.
2. Insert the barbed end of the connector securely into the perfusion tubing.
3. Insert the male luer end of the connector into the desired female luer with a slight twisting motion to insure a secure tight seal.
4. De-bubble the line if it is intended for perfusion.



NovoSci[™]

2021 Airport Rd.
Conroe, TX 77301
Toll Free 1-800-854-0567
Direct 1-281-363-4949

424006-000 Rev. D