

DISPOSABLE PARALLEL “Y” CONNECTOR

R_X Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Intended Use/Indications:

The Parallel “Y” Connector may be used to connect the two caval cannulae to the venous return line and/or make branch connections in the arterial line during cardiopulmonary bypass.

Contraindications:

This Parallel “Y” Connector is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for the cardiopulmonary bypass in general.



Precautions:

Caution: Always de-bubble the Parallel “Y” Connector prior to making patient connection.

Caution: The Parallel “Y” Connector is disposable. Do not reuse, reprocess or re-sterilize.

Direction for Use:

1. Select the Parallel “Y” Connector with the appropriate branch connector sizes for the intended use.
2. Aseptically remove the Parallel “Y” Connector from the peel pouch and transfer it to the sterile field.
3. Push the appropriate connector branches into the appropriate tubes and cannulae making sure all barbs are covered.
4. Completely de-bubble the connector.



NovoSci™

2021 Airport Rd.
Conroe, TX 77301
Toll Free 1-800-854-0567
Direct 1-281-363-4949

424011-000 Rev. F

DISPOSABLE PARALLEL “Y” CONNECTOR

R_X Only

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Intended Use/Indications:

The Parallel “Y” Connector may be used to connect the two caval cannulae to the venous return line and/or make branch connections in the arterial line during cardiopulmonary bypass.

Contraindications:

This Parallel “Y” Connector is not, in itself, a medical treatment. There are, therefore, no absolute contraindications other than those for the cardiopulmonary bypass in general.



Precautions:

Caution: Always de-bubble the Parallel “Y” Connector prior to making patient connection.

Caution: The Parallel “Y” Connector is disposable. Do not reuse, reprocess or re-sterilize.

Direction for Use:

1. Select the Parallel “Y” Connector with the appropriate branch connector sizes for the intended use.
2. Aseptically remove the Parallel “Y” Connector from the peel pouch and transfer it to the sterile field.
3. Push the appropriate connector branches into the appropriate tubes and cannulae making sure all barbs are covered.
4. Completely de-bubble the connector.



NovoSci™

2021 Airport Rd.
Conroe, TX 77301
Toll Free 1-800-854-0567
Direct 1-281-363-4949

424011-000 Rev. F